



# Ryberg Omnia<sup>®</sup>

360 UV-C disinfection

Whitepaper:

Maximizing UV-C Efficiency: A Strategic Approach to Infection Prevention in Hospitals

Autonomous.

Reliable.

Effective.

RYBERG OMNIA: ENGINEERED TO REDEFINE INFECTION PREVENTION

Ryberg is revolutionizing healthcare with cutting-edge UV-C technology. Designed for efficiency, built for impact. Omnia UV-C Disinfection Robot is an autonomous data driven Robot eliminating pathogens with unrivaled precision.

Founded in 2020, Ryberg leads the fight against healthcare-associated infections and antimicrobial resistance. Advanced engineering meets smart automation.



Executive Summary: Hospital-acquired infections (HAIs) pose a significant challenge in healthcare settings, leading to increased patient morbidity, prolonged hospital stays, and higher healthcare costs. Effective disinfection strategies are crucial to mitigating these risks. UV-C light has been widely recognized as a potent method for reducing pathogens on surfaces and in the air. However, the efficiency of UV-C disinfection depends heavily on factors such as distribution, intensity, and exposure duration.

This white paper explores the impact of optimal UV-C distribution on infection prevention in hospitals, focusing on the innovative design features of the Omnia UV-C disinfection robot. Key attributes include 360-degree UV-C distribution, dynamic mobility to reduce shadowing, close-proximity disinfection for improved efficacy, and emission of UV-C wavelengths near the germicidal optimum. Additionally, the integration of a forced airflow system is examined for its role in minimizing airborne pathogen transmission.

## 1. Introduction: The Role of UV-C in Hospital Infection Control

Hospitals are high-risk environments where maintaining a sterile atmosphere is crucial to patient safety. HAIs account for thousands of preventable deaths annually and impose a considerable financial burden on healthcare systems worldwide. Traditional cleaning methods, such as manual wiping with disinfectants, can leave residual pathogens that contribute to infection spread.

UV-C light, with wavelengths between 100 and 280 nm, has been demonstrated to deactivate a broad spectrum of pathogens, including bacteria, viruses, and fungi, by disrupting their DNA and RNA structures. Studies confirm that UV-C technology is highly effective in reducing microbial loads on surfaces and in the air, significantly improving hospital hygiene and patient safety ([Anderson et al., 2017](#); [Rutala & Weber, 2016](#)).

However, the effectiveness of UV-C disinfection is influenced by key factors:

- Intensity: Sufficient UV-C energy is required to ensure microbial inactivation.
- Distribution: Uniform coverage is essential to prevent shadowed areas where pathogens may survive.
- Exposure Duration: Longer exposure times improve efficacy but must be balanced with operational efficiency.

## 2. Design Features of Ryberg's Omnia UV-C Disinfection Robot

Ryberg's Omnia UV-C disinfection robot is engineered to maximize UV-C Distribution, ensuring comprehensive pathogen inactivation while minimizing operational disruptions.

### 2.1 Even UV-C Distribution in All Directions

Key Features:

- Eight Symmetrically Positioned Lamps: A 360° lamp configuration ensures uniform UV-C Distribution, reducing the likelihood of untreated shadowed areas. The robot's 1.2-meter lamps provide extensive high- and low-level coverage, targeting frequently touched surfaces such as bed rails, doorknobs, and medical equipment.





- **Optimized Height:** The robot's height is designed to maximize exposure to critical areas within hospital rooms.

Impact on Infection Control:

- **Comprehensive Coverage:** Ensures high-touch surfaces are thoroughly disinfected.
- **Reduced Shadow Formation:** The symmetric lamp configuration minimizes areas where pathogens could evade exposure.

## 2.2 Dynamic Mobility for Enhanced Disinfection

Key Features:

- **Automated Movement:** Unlike stationary UV-C devices that require manual repositioning, the Omnia robot autonomously navigates through hospital environments.
- **Intelligent Path Planning:** The robot identifies obstacles and optimizes its route to ensure thorough exposure of all surfaces.

Impact on Infection Control:

- **Reduced Shadowing:** Movement allows for exposure of previously blocked surfaces.
- **Labor Efficiency:** Eliminates the need for manual intervention, reducing the risk of user error.

## 2.3 Close-Proximity Disinfection for Improved Efficacy

Key Features:

- **Minimal Distance Operation:** The Omnia robot can function within 15 cm of objects, enhancing UV-C intensity.
- **Narrow Pathway Navigation:** Capable of maneuvering through hospital corridors as narrow as 80 cm.

Impact on Infection Control:

- **Higher Pathogen Inactivation Rates:** Close-range operation increases UV-C exposure, leading to more effective microbial elimination.
- **Energy Efficiency:** Proximity reduces the need for excessive UV-C power, optimizing battery life.

## 3. Optimal Germicidal Wavelength for Maximum Effectiveness

The germicidal efficacy of UV-C is highly dependent on wavelength. Research indicates that 260 nm is the most effective for microbial inactivation, as it corresponds to peak DNA absorption ([Dai et al., 2012](#)).

Key Features:

- **Emission at 253.7 nm:** The Omnia robot operates at a wavelength close to the germicidal optimum.



- Non-Penetration of Glass: Unlike broad-spectrum pulsed xenon UV light, 253.7 nm UV-C does not pass through glass, ensuring focused disinfection.

Impact on Infection Control:

- Maximized Pathogen Elimination: Ensures high disinfection efficiency.
- Enhanced Safety: Prevents unintended UV-C exposure in adjacent areas.

#### 4. Forced Airflow for Enhanced Pathogen Control

Airborne pathogens present a significant infection risk, particularly in enclosed hospital settings. The Omnia UV-C robot incorporates a forced airflow system that enhances airborne pathogen control.

Key Features:

- Upward and Downward Airflow Circulation: Channels air through UV-C exposure zones before recirculating it into the room.
- Minimized Airborne Spread: Designed to limit the transmission of pathogens through air currents.

Impact on Infection Control:

- Reduction of Airborne Pathogens: Improves air quality and lowers transmission risks.
- Comprehensive Disinfection: Enhances the effectiveness of surface and air pathogen removal.

#### 5. Conclusion: The Future of UV-C Disinfection in Hospitals

Ryber's Omnia UV-C disinfection robot represents a significant advancement in infection prevention. Its innovative design—featuring optimized UV-C distribution, autonomous mobility, close-proximity disinfection, and forced airflow—ensures superior efficacy in pathogen elimination.

As hospitals continue to confront infection control challenges, the integration of advanced disinfection technologies like Omnia UV-C robot will play a pivotal role in reducing HAIs, improving patient outcomes, and enhancing overall hospital safety. Ongoing research and technological advancements in UV-C disinfection will further refine and expand its applications, ensuring that healthcare facilities remain at the forefront of infection prevention.

#### References

- Anderson, D. J., et al. (2017). "Enhanced Terminal Room Disinfection and Acquisition of Multidrug-Resistant Organisms and *Clostridium difficile* in Healthcare Settings." JAMA. [<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC5626443/>]
- Rutala, W. A., & Weber, D. J. (2016). "CDC Guidelines for Disinfection and Sterilization in Healthcare Facilities." [<https://www.cdc.gov/infectioncontrol/pdf/guidelines/disinfection-guidelines-H.pdf>]
- Dai, T., et al. (2012). "Ultraviolet C Irradiation: An Alternative Antimicrobial Approach to Localized Infections?" [<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC3486633/>]

